


Easier & Safer Usage

Combined Torque Wrench



This is a **brand new type of Combined torque wrench** that finishes screwing and disassembling any fixtures and abutments with a single device in a one-step process.

- **It does not require an extra hand driver**



Convenience

Compact and light-weight, a wide field of view is possible even with a limited mouth opening.




Safety

Completely safe operation with the innovative structure of an all-in-one driver & torque wrench.

- **Specifications**

Torque Scale	20~50 Ncm
Cat No.	Product Description
DFCTW-W-N5012-RH	Denflex Combined Torque Wrench without Driver (Only Torque Wrench)
DFCTW-D-H1213-OS	1.2mm Hex Driver for Denflex Combined Torque Wrench, Length 13mm
DFCTW-D-H1213-DT	1.25mm Hex Driver for Denflex Combined Torque Wrench, Length 13mm
DFCTW-D-R1616-ST	SCS Driver for Denflex Combined Torque Wrench, Length 16mm
DFCTW-D-R1516-NB	Unigrip Driver for Denflex Combined Torque Wrench, Length 16mm
DFFAA-D-H1217-OS	1.2mm Hex Driver, Length 17mm
DFFAA-D-H1217-DT	1.25mm Hex Driver, Length 17mm
DFFAA-D-H2618-OS	Ball Driver, 2.6mm Hex, Length 18mm

WHY Denflex Combined Torque Wrench?



• What is a torque wrench in general?



The initial use of a hand driver




The additional uses of a torque wrench

The existing torque wrenches require a separate hand driver, and it is easy to lose the hand driver while operating in the mouth, possibly causing the patient to swallow it. Patients with a limited mouth opening or surgeons with large hands are vulnerable to great discomfort especially during surgery on posterior teeth.

• Features of the Denflex Combined Torque Wrench

An all-in-one driver & torque wrench




No extra uses of a hand driver.



A single torque wrench from the initial tightening through to the finish.



A wide field of view is guaranteed even with a limited mouth opening.



It has a compact and lightweight structure.

Denflex
www.denflex.kr

Denflex Co., Ltd.
#708, Byucksan Digital Valley I, 24, Digital-ro 27-gil,
Guro-gu, Seoul, Republic of Korea
Tel. 02-866-0617 **Fax.** 02-866-0618 **Email.** info@denflex.kr

